


Unter Hafßen 20 · D-78532 Tuttlingen
Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50
e-mail: sales@rz-medizintechnik.com
<http://www.rz-medizintechnik.com>




UROLOGY

Cystoscopy



UROLOGY CYSTOSCOPY



Flexible Instruments

Biopsy Forceps, Graspers and Scissors are available in 7 Charr. and 9 Charr. with 400mm working length



Cysto-Urethroscope Sheaths

are colour coded and available in 17, 19, 20, 22 and 25 Charr. All sheaths come with atraumatic obturators and 2 stop cocks for suction & irrigation.



Bridges

Are available with 2, with 1 or without working channels. For use with all Cystoscope Sheaths and flexible instruments.



Telescopes

Come in 4.0mm and 2.9mm diameter and 0°, 12°, 30° and 70° angle of view. Distal and Proximal saphire glasses protect the optical system from damages.



Rigid Optical Forceps

With anatomical designed handle and a wide range of different jaws, e.g. biopsy forceps, grasper, punch and scissor



Albarran Deflecting Mechanisms

Are used to deflect flexible instruments like electrodes or biopsy forceps to reach anatomically difficult areas. The built-in ratchet mechanism keeps the instrument in place during surgery.




HF Electrodes

In 5, 7 and 9 Charr. are available in ball, needle and loop style for all kinds of Coagulation. Suitable cables for all generators are available.

Your Advantages of the RZ Medizintechnik GmbH Cystoscope Set are:

- Colour coding for easy identification
- Wide angle telescopes for a maximum field of view
- Only Stainless Steel and medical grade plastic components
- Fully autoclavable
- Unique and lightweight designs for maximum comfort
- Compatible to other brands
- Repair & Repair-Exchange Service for all major telescope and instrument brands

BASIC CYSTOSCOPE SET



| | |
|-------------|---|
| 253-804-030 | Cystoscope Telescope, 4mm, 30° |
| 253-300-017 | Cystoscope Sheath 17 Charr. incl. Obturator |
| 253-300-022 | Cystoscope Sheath 22 Charr. incl. Obturator |
| 253-000-102 | Bridge with 2 channel |
| 253-000-105 | Albarran deflecting mechanism with 2 channel |
| 253-420-100 | Optical Biopsy Forceps, double action, round cup |
| 253-400-307 | Flexible Grasping Forceps, 7 Charr. |
| 253-400-407 | Flexible Biopsy Forceps, 7 Charr. |
| 253-370-005 | Flexible Ball Electrode, 5 Charr. |
| 253-610-004 | HF Cable with 4mm plug |
| 600-464-229 | Plastic Container with silicone mat and transparent lid |