


Unter Haßlen 20 · D-78532 Tuttlingen
Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50
e-mail: sales@rz-medizintechnik.com
http://www.rz-medizintechnik.com



HYSTEROSCOPY

Continuous High-Flow Sheath




03/10

PERFECT CONTINUOUS FLOWRATES WITH THE NEW DESIGNED RZ HYSTEROSCOPY SHEATH

This single-piece sheath for use with 2.9mm telescopes offers advanced flowrates, even with inserted instrument.

| 351-400-055

MIDI Hysteroscopic Continuous Flow sheath
5.5mm, 5 Charr. Instrument channel, for
2.9mm, 30° scope




HD
ready

| 351-829-030 Hysteroscope, 30° direction of view, 2,9mm diameter

5 Charr. Hysteroscopy Instruments



| 351-100-305 Scissors, blunt tips, single action



| 351-118-305 Scissors, sharp tips, single action



| 351-116-305 Grasping Forceps with cup, serrated, double action



| 351-120-305 Biopsy Forceps, double action



| 351-126-305 Biopsy Punch Forceps, single action



| 351-124-105 Myom fixation screw, flexible, hexagonal handle



| 351-370-105 Bipolar Button Electrode, 36cm, 5 Charr. , flexible



| 351-360-005 Bipolar Needle Electrode, 36cm, 5 Charr. , flexible



| 700-350-030 Bipolar Connecting Cable, 3m length

| 700-350-050 Bipolar Connecting Cable, 5m length

other cables on request